



PRODUCT SPECIFICATIONS

Date: 2011-08-12

PRODUCT NAME : Anti-human MPO 1702 SPTN-5

PRODUCT SPECIFICITY : Antibody recognizes human Myeloperoxidase

PRODUCT CODE : 100267

PRODUCT BUFFER : 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.1 % NaN₃ as a preservative

SHELF LIFE AND STORAGE : 12 months at 2-8 °C

ANALYTE DESCRIPTION : Myeloperoxidase (MPO) is a peroxidase enzyme (EC 1.11.1.7) most abundantly

present in neutrophil granulocytes. The 150-kDa MPO is a dimer consisting of two 15-kDa light chains and two variable-weight glycosylated heavy chains bound to a prosthetic heme group. Myeloperoxidase (MPO) is an enzymatic mediator of several inflammatory cascades and higher serum levels have been associated with increased risk of adverse cardiovascular events. As a result, MPO and its downstream inflammatory pathways represent attractive targets for both prognostic and therapeutic

intervention in the prophylaxis of atherosclerotic cardiovascular disease.

PARAMETERS TESTED FROM EACH LOT

PRODUCT APPEARANCE : Clear liquid, may turn opaque during storage

PRODUCT CONCENTRATION : 5.0 mg/ml (+/- 10 %)

PRODUCT ACTIVITY : 80 – 120 % compared to reference in an IFMA-test

IEF-RANGE : 6.5 – 7.4

PURITY : ≥ 95 %

PARAMETERS DETERMINED ONLY DURING PRODUCT R&D PHASE

CLASS AND SUBCLASS : IgG1

ASSOCIATION CONSTANT : 8 x 10⁵ 1/Ms

DISSOCIATION CONSTANT : 2 x 10⁻³ 1/s

AFFINITY CONSTANT : $K_D = 2 \times 10^{-9} \text{ M}; K_A = 5 \times 10^8 \text{ 1/M}$

DETERMINATION METHOD : SPR analysis (ProteOn XPR36)

ANTIGEN : Myeloperoxidase (Lee Biosolutions, cat.no.426-10, lot.no. L44340)

CROSS-REACTIVITIES : N/D



CARDIOVASCULAR

medix_{MAB}^m

EPITOPE : N/D

EPITOPE GROUP : -

Two antibodies binding to the same, or closely located epitopes, belong to the same group and hence cannot be used as a pair in a sandwich assay. Epitope group numbering does not give any detailed information where the epitope is located.

PAIR RECOMMENDATIONS : SOLID LABEL

1702 1703

Please note that pair recommendations are based on results obtained in our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations should be taken only as a directive.

PRODUCT STABILITY : TEMPERATURE, DAYS RESULT

Please note that the shelf life given on page one is based on real time stability testing at +2-8 °C in the product buffer.

-70 °C, 21 days

-20 °C, 21 days

+4 °C, 21 days

OK

+25 °C, 21 days

OK

+35 °C, 21 days

OK

-45 °C, 7 days

OK

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. The maximum duration of the test is 21 days, except for the +45 $^{\circ}\text{C}$ only 7 days.

pH, 14 DAYS, +4°C	RESULT
5.0	OK
6.0	OK
7.0	OK
8.0	ОК

Stability testing is performed to see whether pH affects the antigen binding, charge or composition of the antibody during 14 days at +4 $^{\circ}$ C.

MISCELLANEOUS : -

REFERENCES : -

Legal disclaimer

MedixMAB monoclonal antibodies meet their specifications if transported, stored and used according to the instructions for use. See also separate Notes on use and storage of monoclonal antibodies. If not otherwise agreed in writing all products are sold under General Delivery Terms. Other terms may apply if MedixMAB monoclonal antibodies are purchased from an official distributor of Medix Biochemica products. MedixMAB is a registered trademark of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission. Further information on products and methods is available from medixMAB website and extranet.