



PRODUCT SPECIFICATIONS

Date: 2012-03-22

PRODUCT NAME	:	Anti-h LH 5302 SP-1
PRODUCT SPECIFICITY	:	Antibody recognizes human Luteinizing Hormone (Lutropin), and its beta-subunit
PRODUCT CODE	:	100018
PRODUCT BUFFER	:	0.9 % NaCl, 0.1 % NaN₃ as a preservative
SHELF LIFE AND STORAGE	:	24 months from manufacturing at 2-8 °C
ANALYTE DESCRIPTION	:	In both males and females, LH is essential for reproduction. In females FSH initiates follicular growth and at the time of the maturation of the follicle the estrogen rise leads to a release of LH over a 24-48 hour period. This 'LH surge' triggers ovulation thereby not only releasing the egg, but also initiating the conversion of the residual follicle into a corpus luteum that, in turn, produces progesterone to prepare the endometrium for a possible implantation. LH is necessary to maintain luteal function for the first two weeks. In case of a pregnancy luteal function will be further maintained by the action of hCG from the newly established pregnancy. In the male, LH acts upon the Leydig cells of the testis and is responsible for the production of

PARAMETERS TESTED FROM EACH LOT

PRODUCT APPEARANCE	:	Clear liquid
PRODUCT CONCENTRATION	:	1.00 mg/ml (+/- 10 %)
PRODUCT ACTIVITY	:	80-120 % compared to reference in an IFMA-test
IEF-RANGE	:	5.1 - 5.7
PURITY	:	≥ 95 %

testosterone.

PARAMETERS DETERMINED ONLY DURING PRODUCT R&D PHASE

CLASS AND SUBCLASS	:	lgG1		
ASSOCIATION CONSTANT	:	1.3 x 10 ⁶ 1/Ms		
DISSOCIATION CONSTANT	:	2.6 x 10 ⁻⁵ 1/s		
AFFINITY CONSTANT	:	KA = 4.8 x 10 ¹⁰ 1/M; KD = 2.1 x 10 ⁻¹¹ M (= 21 pM)		
DETERMINATION METHOD	:	SPR analysis (ProteOn XPR36)		
ANTIGEN	:	LH, Scripps Laboratories (Cat No L0815 Lot#2360102)		
CROSS-REACTIVITIES	:	$ \begin{array}{llllllllllllllllllllllllllllllllllll$		

TSH 0.02 % (Scripps Laboratories Cat No T0114 Lot 181711)

Oy Medix Biochemica Ab Asematie 13 FI-02700 Kauniainen Finland Telephone +358-9-547 680 Fax +358-9-505 3441 Email medix@medixbiochemica.com Public website www.medixbiochemica.com Extranet

www.medixbiochemica.com/extra Request your personal login ID by emailing to medix@medixbiochemica.com



FERTILITY

EPITOPE	:			
EPITOPE GROUP	:	Beta 1 as showed in Nilsson et al. (2001)		
		Two antibodies binding to the same, or closely located epitopes, belong to the same group and hence cannot be used as a pair in a sandwich assay. Epitope group numbering does not give any detailed information where the epitope is located.		
PAIR RECOMMENDATIONS	:	SOLID	LABEL	
		5302	5301	
		5301	5302	
		Please note that pair recommendations are based on results obtained in our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations should be taken only as a directive.		
PRODUCT STABILITY	:	TEMPERATURE, DAYS	RESULT	
		-70 °C, 21 days	N/D	
		-20 °C, 21 days	N/D	
Please note that the shelf life given on		+4 °C, 21 days	N/D	
page one is based on real time stability		+25 °C, 21 days	N/D	
testing at +2-8 °C in the product buffer.		+35 °C, 7 days	N/D	
		+35 °C, 21 days	N/D	
		+45 °C, 3 days	N/D	
		+45 °C , 7 days	N/D	
		Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. The maximum duration of the test is 21 days, except for the +45 °C only 7 days.		
		pH, 14 DAYS, +4°C	RESULT	
		5.0	N/D	
		6.0	N/D	
		7.0	N/D	
		8.0	N/D	
		Stability testing is performed to see whether composition of the antibody during 14 days	er pH affects the antigen binding, charge or s at +4 $^{\circ}$ C.	
MISCELLANEOUS	:	In Nilsson et al. (2001) authors performed an epitope mapping analysis of 30 different monoclonal antibodies for LH. In their analysis it was shown that 5302 binds to an epitope named as beta 1. Antibodies reacting with beta 1 epitope did not react with FSH. They also showed that 5302 reacts with intact LH as well as with the beta-subunit of LH while having some cross-reactivity with hCG but no cross-reactivity with TSH or FSH. Clone 5302 reacted also with the variant-LH.		
REFERENCES		Federici, M.M., Fraser, R., Lundqvist, C., and Lankford, J.C., (1982) Production and characterization of monoclonal antibodies to human luteinizing hormones. Fed. Proc. 41		
		Nilsson, C., Seppälä, M., and Pettersson, of human luteinizing hormone with special J.Endocrinol. 168:10-116		

Legal disclaimer

MedixMAB monoclonal antibodies meet their specifications if transported, stored and used according to the instructions for use. See also separate *Notes on use and storage of monoclonal antibodies*. If not otherwise agreed in writing all products are sold under *General Delivery Terms*. Other terms may apply if MedixMAB monoclonal antibodies are purchased from an official distributor of Medix Biochemica products. MedixMAB is a registered trademark of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission. Further information on products and methods is available from medixMAB website and extranet.

Oy Medix Biochemica Ab Asematie 13 FI-02700 Kauniainen Finland **Telephone** +358-9-547 680 **Fax** +358-9-505 3441

Email medix@medixbiochemica.com Public website www.medixbiochemica.com

Extranet

www.medixbiochemica.com/extra Request your personal login ID by emailing to medix@medixbiochemica.com