


**INFECTIOUS
DISEASES**
PRODUCT SPECIFICATIONS

Date: 2010-08-03

PRODUCT NAME	: Anti-Chlamydia 6701 SP-1
PRODUCT SPECIFICITY	: Antibody recognizes <i>Chlamydia trachomatis</i> LPS
PRODUCT CODE	: 100068
PRODUCT BUFFER	: 0.9 % NaCl, 0.1 % NaN ₃ as a preservative
SHELF LIFE AND STORAGE	: 24 months from manufacturing at 2-8 °C
ANALYTE DESCRIPTION	: <i>C. trachomatis</i> is a obligate intracellular pathogen and can cause numerous disease states in both men and women. Both sexes can display urethritis, proctitis, trachoma, and infertility. The bacterium can cause prostatitis and epididymitis in men. In women, cervicitis, pelvic inflammatory disease (PID), ectopic pregnancy, and acute or chronic pelvic pain are frequent complications. <i>C. trachomatis</i> is also an important neonatal pathogen, where it can lead to infections of the eye (trachoma) and pulmonary complications.

PARAMETERS TESTED FROM EACH LOT

PRODUCT APPEARANCE	: Clear liquid
PRODUCT CONCENTRATION	: 1.00 mg/ml (+/- 10 %)
PRODUCT ACTIVITY	: 80-120 % compared to reference in an IFMA-test
IEF-RANGE	: 5.6 - 6.7
PURITY	: ≥ 95 %

PARAMETERS DETERMINED ONLY DURING PRODUCT R&D PHASE

CLASS AND SUBCLASS	: IgG ₁
ASSOCIATION CONSTANT	:
DISSOCIATION CONSTANT	:
AFFINITY CONSTANT	: n.d.
DETERMINATION METHOD	:
ANTIGEN	: For immunization Chlamydia LGV2 strain
CROSS-REACTIVITIES	: <i>C. pneumoniae</i> <i>C. psittaci</i> <i>S. aureus</i> (unspecific binding to protein A)

Cross-reactivity with the following micro-organisms below 5%

A. lwoffii, *A. baumannii*, *B. fragilis*, *C. albicans*, *C. freundii*, *C. xerosis*, *E. aerogenes*, *E. coli*, *Fusobacterium*, *G. vaginalis*, *K. pneumonia*, *L. acidophilus*, *M. hominis*, *M. smegmatis*, *N. gonorrhoeae*, *Peptococcus*, *Peptostreptococcus*, *Propionebacterium*, *P. vulgaris*, *R. aeruginosa*, *S. cerevisiae*, *S. marcescens*, *S. faecalis*, *S. pyogenes*, *T. globrata*, *U. urealyticum*, *S. minnesota*

EPITOPE	:	Chlamydial LPS KDO-trisaccharide																												
EPITOPE GROUP	:	Two antibodies binding to the same, or closely located epitopes, belong to the same group and hence cannot be used as a pair in a sandwich assay. Epitope group numbering does not give any detailed information where the epitope is located.																												
PAIR RECOMMENDATIONS	:	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left;">SOLID</th> <th style="text-align: left;">LABEL</th> </tr> </thead> <tbody> <tr> <td>6709</td> <td>6701</td> </tr> </tbody> </table> <p>Please note that pair recommendations are based on results obtained in our laboratory and equally good results can be obtained using also other pairs and therefore these recommendations should be taken only as directive.</p>	SOLID	LABEL	6709	6701																								
SOLID	LABEL																													
6709	6701																													
PRODUCT STABILITY	:	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left;">TEMPERATURE, DAYS</th> <th style="text-align: left;">RESULT</th> </tr> </thead> <tbody> <tr><td>-70 °C, 21 days</td><td></td></tr> <tr><td>-20 °C, 21 days</td><td></td></tr> <tr><td>+4 °C, 21 days</td><td></td></tr> <tr><td>+25 °C, 21 days</td><td></td></tr> <tr><td>+35 °C, 7 days</td><td></td></tr> <tr><td>+35 °C, 21 days</td><td></td></tr> <tr><td>+45 °C, 3 days</td><td></td></tr> <tr><td>+45 °C, 7 days</td><td></td></tr> </tbody> </table> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left;">pH, 14 DAYS, +4 °C</th> <th style="text-align: left;">RESULT</th> </tr> </thead> <tbody> <tr><td>5.0</td><td></td></tr> <tr><td>6.0</td><td></td></tr> <tr><td>7.0</td><td></td></tr> <tr><td>8.0</td><td></td></tr> </tbody> </table>	TEMPERATURE, DAYS	RESULT	-70 °C, 21 days		-20 °C, 21 days		+4 °C, 21 days		+25 °C, 21 days		+35 °C, 7 days		+35 °C, 21 days		+45 °C, 3 days		+45 °C, 7 days		pH, 14 DAYS, +4 °C	RESULT	5.0		6.0		7.0		8.0	
TEMPERATURE, DAYS	RESULT																													
-70 °C, 21 days																														
-20 °C, 21 days																														
+4 °C, 21 days																														
+25 °C, 21 days																														
+35 °C, 7 days																														
+35 °C, 21 days																														
+45 °C, 3 days																														
+45 °C, 7 days																														
pH, 14 DAYS, +4 °C	RESULT																													
5.0																														
6.0																														
7.0																														
8.0																														
MISCELLANEOUS	:	Clone 6701 binds to chlamydial lipopolysaccharide (LPS) and therefore recognizes all chlamydias.																												
REFERENCES	:	Tirola, T., Jakkola, A., Bloigu, A., Paldanius, M., Sinisalo, J., Nieminen, M.S., Silvennoinen-Kassinen, S., Saikku, P., Jauhainen, M. and Leinonen, M. (2006) Novel enzyme immunoassay utilizing lipopolysaccharide-binding protein as a capture molecule for the measurement of chlamydial lipopolysaccharide in serum. <i>Diag. Microbiol. Inf. Dis.</i> 54:7-12																												

Legal disclaimer

MedixMAB monoclonal antibodies meet their specifications if transported, stored and used according to the instructions for use. See also separate *Notes on use and storage of monoclonal antibodies*. If not otherwise agreed in writing all products are sold under *General Delivery Terms*. Other terms may apply if MedixMAB monoclonal antibodies are purchased from an official distributor of Medix Biochemica products. MedixMAB is a registered trademark of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission. Further information on products and methods is available from medixMAB website and extranet.