


CARDIAC
PRODUCT SPECIFICATIONS
Date: 2010-05-24

| | | |
|-------------------------------|---|---|
| PRODUCT NAME | : | Anti-h cTnI 9703 SPTN-5 |
| PRODUCT SPECIFICITY | : | Antibody recognizes human Cardiac Troponin I |
| PRODUCT CODE | : | 100181 |
| PRODUCT BUFFER | : | 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.1 % NaN ₃ as a preservative |
| SHELF LIFE AND STORAGE | : | 6 months from manufacturing at 2-8 °C |
| ANALYTE DESCRIPTION | : | Troponins form a complex (Troponin C, I and T) of three regulatory proteins that are integral to muscle contraction in skeletal and cardiac muscle, but not smooth muscle. Troponin levels can be used as a test of several different heart disorders, including myocardial infarction. The measurement of serum cTnI and cTnT is superior in terms of sensitivity and specificity to cardiac muscle enzyme measurements in the identification of cardiac muscle damage. Raised cardiac troponin concentrations are now accepted as the standard biochemical marker for the diagnosis of myocardial infarction. |

PARAMETERS TESTED FROM EACH LOT

| | | |
|------------------------------|---|--|
| PRODUCT APPEARANCE | : | Clear liquid |
| PRODUCT CONCENTRATION | : | 5.0 mg/ml (+/- 10 %) |
| PRODUCT ACTIVITY | : | 80-120 % compared to reference in an IFMA-test |
| IEF-RANGE | : | 5.2 - 6.7 |
| PURITY | : | ≥ 95 % |

PARAMETERS DETERMINED ONLY DURING PRODUCT R&D PHASE

| | | |
|------------------------------|---|--|
| CLASS AND SUBCLASS | : | IgG ₁ |
| ASSOCIATION CONSTANT | : | |
| DISSOCIATION CONSTANT | : | |
| AFFINITY CONSTANT | : | 1.0 x 10 ⁹ l/mol |
| DETERMINATION METHOD | : | Biacore |
| ANTIGEN | : | Troponin complex, HyTest (Cat No 8T53 Lot 06/07) |
| CROSS-REACTIVITIES | : | |

EPITOPE : Between residues 27-73

EPITOPE GROUP :

Two antibodies binding to the same, or closely located epitopes, belong to the same group and hence cannot be used as a pair in a sandwich assay. Epitope group numbering does not give any detailed information where the epitope is located.

PAIR RECOMMENDATIONS :

| SOLID | LABEL |
|-------|-------|
| 9707 | 9703 |

Please note that pair recommendations are based on results obtained in our laboratory and equally good results can be obtained using also other pairs and therefore these recommendations should be taken only as directive.

PRODUCT STABILITY :

| TEMPERATURE, DAYS | RESULT |
|--------------------|--------|
| -70 °C, 21 days | |
| -20 °C, 21 days | |
| +4 °C, 21 days | |
| +25 °C, 21 days | |
| +35 °C, 7 days | |
| +35 °C, 21 days | |
| +45 °C, 3 days | |
| +45 °C, 7 days | |
| pH, 14 DAYS, +4 °C | RESULT |
| 5.0 | |
| 6.0 | |
| 7.0 | |
| 8.0 | |

MISCELLANEOUS :

REFERENCES :

Legal disclaimer

MedixMAB monoclonal antibodies meet their specifications if transported, stored and used according to the instructions for use. See also separate *Notes on use and storage of monoclonal antibodies*. If not otherwise agreed in writing all products are sold under *General Delivery Terms*. Other terms may apply if MedixMAB monoclonal antibodies are purchased from an official distributor of Medix Biochemica products. MedixMAB is a registered trademark of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission. Further information on products and methods is available from medixMAB website and extranet.