



PRODUCT SPECIFICATION

CATALOG NUMBER: J87500015 = 2-001-A

PRODUCT DESCRIPTION:	Mouse IgG
SOURCE:	Normal mouse serum
APPEARANCE:	Clear, precipitate-free liquid
PURITY:	>97% by SDS-PAGE
BUFFER:	10mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS)
PRESERVATIVE:	0.05% (w/v) sodium azide
MICROBIAL:	0 CFU/mL
STORAGE TEMPERATURE:	2-8°C
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative. Sodium azide may react with lead plumbing to form explosive metal azides. Flush with large amounts of water during disposal. There is no known biohazard associated with this product.