



## PRODUCT SPECIFICATION

**CATALOG NUMBER: 2-025-C**

PRODUCT:	Purified 4C5 Cardiac Troponin T Monoclonal Antibody
IMMUNOGEN:	Highly purified human cardiac troponin T subunit
FORM:	Liquid
PURITY:	>95% by SDS-PAGE. Purified by protein A agarose chromatography.
BUFFER:	10 mM Sodium Phosphate, 0.85% (w/v) Sodium Chloride, pH 7.2
PRESERVATIVE:	0.05% (w/v) sodium azide
STORAGE:	2-8 degrees C