

PRODUCT SPECIFICATIONS

Date: 2015-05-12

PRODUCT NAME	:	Anti-h IL-6 2703 SPTN-5
PRODUCT SPECIFICITY	:	Human Interleukin-6 (IL-6)
PRODUCT CODE	:	100328
PRODUCT BUFFER	:	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative.
SHELF LIFE AND STORAGE	:	24 months at 2-8 °C
ANALYTE DESCRIPTION	:	IL-6 is an interleukin that acts as both a pro-inflammatory and anti-inflammatory cytokine. It is secreted by T cells and macrophages to stimulate immune response to trauma, especially burns or other tissue damage leading to inflammation. IL-6 is relevant to many disease processes such as diabetes, atherosclerosis, depression, Alzheimer's Disease, systemic lupus erythematosus, cancer, and rheumatoid arthritis. There is an interest in developing anti-IL-6 agents as therapy against many of these diseases. The first such is tocilizumab which has been approved for rheumatoid arthritis.

PARAMETERS TESTED FROM EACH LOT

PRODUCT APPEARANCE	:	Clear liquid, may turn opaque during storage
PRODUCT CONCENTRATION	:	5.0 mg/ml (+/- 10 %)
PRODUCT ACTIVITY	:	80 – 120 % compared to reference in an IFMA-test
IEF-RANGE	:	6.4 – 7.3
PURITY	:	≥ 95 %

PARAMETERS DETERMINED ONLY DURING PRODUCT R&D PHASE

CLASS AND SUBCLASS	:	IgG ₁
ASSOCIATION CONSTANT	:	1.0 x 10 ⁵ 1/Ms
DISSOCIATION CONSTANT	:	5.6 x 10 ⁻⁵ 1/s
AFFINITY CONSTANT	:	KA = 1.8 x 10 ⁹ 1/M; KD = 5.6 x 10 ⁻¹⁰ M (=0.56 nM)
DETERMINATION METHOD	:	SPR analysis (BioRad ProteOn XPR36)
ANTIGEN	:	IL-6 (Creative Biomart cat.no. IL6-12H, lot#714033)
CROSS-REACTIVITIES	:	N/D

EPITOPE : N/D

EPITOPE GROUP :

Two antibodies binding to the same, or closely located epitopes, belong to the same group and hence cannot be used as a pair in a sandwich assay. Epitope group numbering does not give any detailed information where the epitope is located.

PAIR RECOMMENDATIONS	SOLID	LABEL
	2703	2704
	2704	2703

Please note that pair recommendations are based on results obtained in our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations should be taken only as a directive.

PRODUCT STABILITY	TEMPERATURE, DAYS	RESULT
Please note that the shelf life given on page one is based on real time stability testing at +2-8 °C in the product buffer.	-70 °C, 21 days	OK
	-20 °C, 21 days	OK
	+4 °C, 21 days	OK
	+25 °C, 21 days	OK
	+35 °C, 21 days	OK
	+45 °C, 7 days	OK

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. The maximum duration of the test is 21 days, except for the +45 °C only 7 days.

pH, 14 DAYS, +4 °C	RESULT
5.0	OK
6.0	OK
7.0	OK
8.0	OK

Stability testing is performed to see whether pH affects the antigen binding, charge or composition of the antibody during 14 days at +4 °C.

MISCELLANEOUS : -

REFERENCES : -

Legal disclaimer

MedixMAB monoclonal antibodies meet their specifications if transported, stored and used according to the instructions for use. See also separate *Notes on use and storage of monoclonal antibodies*. If not otherwise agreed in writing all products are sold under *General Delivery Terms*. Other terms may apply if MedixMAB monoclonal antibodies are purchased from an official distributor of Medix Biochemica products. MedixMAB is a registered trademark of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission. Further information on products and methods is available from medixMAB website and extranet.