



## PRODUCT SPECIFICATION

**CATALOG NUMBER: 3-030**

PRODUCT DESCRIPTION:	Purified Rabbit Skeletal Muscle Troponin, C subunit
SOURCE:	Rabbit skeletal muscle
FORM:	Liquid
PURITY:	>95% by SDS-PAGE
BUFFER:	10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS)
PRESERVATIVE:	0.05% (w/v) sodium azide
STORAGE TEMPERATURE:	2-8°C
HAZARD/BIOHAZARD:	Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal. There is no known biohazard associated with this product.