


PRODUCT SPECIFICATIONS

2015-09-01

PRODUCT NAME	: Trypsinogen-2 (Trpg-2)
ANALYTE DESCRIPTION	: Trypsinogen-2 is the precursor form of the pancreatic enzyme trypsin. Determination of urine trypsinogen-2 is a useful test to detect acute pancreatitis and to evaluate disease severity.
PRODUCT CODE	: 610020
PRODUCT DESCRIPTION	: Recombinant human protein, expressed in <i>E. coli</i> UniProt accession # P07478, amino acids 16-247 Predicted molecular weight: 24.9 kD
FORMULATION	: Chromatographically purified recombinant protein, lyophilized from 0.2 µm filtered solution.
QUANTITY	: 100 µg
CONCENTRATION	: 0.5 mg/ml when reconstituted with 200 µl of deionised water
BUFFER	: Tris buffered saline, pH 7.4, containing 3 % sucrose as stabilizer
RECONSTITUTION	: Reconstitute lyophilized protein with 200 µl of deionised water
PURITY	: >90 % by SDS-PAGE



STORAGE	: Store at 2-8 °C
REACTIVITY WITH MedixMABS	: MedixMAB 8603 Anti-h Trypsinogen-2: + MedixMAB 8607 Anti-h Trypsinogen-2: +

Legal disclaimer

MedixMAB monoclonal antibodies, MedixMAB Fab fragments and MedixAntigens meet their specifications if transported, stored and used according to the instructions for use. If not otherwise agreed in writing all products are sold under *General Delivery Terms*. Further information on products is available on the internet.

MedixMAB is a registered trademark of Medix Biochemica.

Oy Medix Biochemica Ab

HEADQUARTERS

Klovinpellontie 3
FI-02180 Espoo, Finland

MANUFACTURING SITE

Noljakantie 13
FI-80130 Joensuu, Finland

medix@medixbiochemica.com

www.medixbiochemica.com

VAT reg. no. FI14631532