



PRODUCT SPECIFICATION

CATALOG NUMBER: A03010

| | |
|---------------------|--|
| PRODUCT: | Monoclonal Antibody to SARS |
| FORM: | 0.02 M phosphate buffer, 0.25M NaCl |
| PRESERVATIVE: | 0.1% sodium azide |
| SUBCLASS: | IgG _{2a} |
| STORAGE CONDITIONS: | 2-8°C |
| SPECIFICITY: | This mab reacts with SARS associated coronavirus (SARS-CoV) S proteins. |
| APPLICATIONS: | Immunofluorescence, ELISA, and Western Blot |
| HAZARD/BIOHAZARD: | This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product. |