



PRODUCT SPECIFICATION

CATALOG NUMBER: A03020

PRODUCT:	Monoclonal Antibody to SARS
FORM:	0.02 M phosphate buffer, 0.25M NaCl
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG _{2a}
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	This mab reacts with SARS associated coronavirus (SARS-CoV) S proteins.
APPLICATIONS:	Immunofluorescence, ELISA, and Western Blot
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.