



PRODUCT SPECIFICATION

catalog number:	A03070
product description:	SARS Coronavirus Monoclonal Antibody
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	0.01 M PBS, pH 7.2
preservative:	0.1% NaN ₃
isotype:	IgG _{2b}
concentration:	[lot specific] mg/mL
specificity:	Specific for nucleoprotein of the SARS virus (IFA, ELISA & WB) AND COVID-19 (ELISA)
notes/applications:	ELISA, IFA and Western Blot
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.