



## PRODUCT SPECIFICATION

catalog number:	<b>A03091</b>
product description:	SARS Coronavirus Monoclonal Antibody
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	0.01 M PBS, pH 7.2
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG
concentration:	[lot specific] mg/mL
specificity:	Specific for nucleoprotein of the SARS virus (IFA, ELISA) and COVID-19 (ELISA)
notes/applications:	ELISA, IFA
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.