



## PRODUCT SPECIFICATION

catalog number:	<b>A03093</b>
product description:	COVID-19 Nucleocapsid Monoclonal Antibody
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	0.01 M PBS, pH 7.2
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
specificity:	Specific for the nucleocapsid of COVID-19 virus.
cross-reactivity:	Cross reacts with SARS. Non-reactive with Influenza A, Influenza B, Adenovirus & MERS.
pairing:	This antibody can be used in sandwich ELISA as the capture antibody in conjunction with A03094 or A03095 or A03096 as the detection antibodies.
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.