



PRODUCT SPECIFICATION

catalog number:	A03094
product description:	COVID-19 Nucleocapsid Monoclonal Antibody
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	0.01 M PBS, pH 7.2
preservative:	0.1% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	Specific for the nucleocapsid of COVID-19 virus.
cross-reactivity:	Cross reacts with SARS. Non-reactive with Influenza A, Influenza B, Adenovirus & MERS.
pairing:	This antibody can be used in sandwich ELISA as the detection antibody in conjunction with A03093 or A03097 as the capture antibodies.
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.