



## PRODUCT SPECIFICATION

catalog number:	<b>A03099</b>
product description:	SARS-CoV-2 Nucleocapsid Monoclonal Antibody
lot number:	[lot specific]
source:	Murine
form:	Liquid
purity:	≥95% by SDS-PAGE
buffer:	PBS pH 7.4
preservative:	0.09% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
specificity:	Specific for SARS-CoV-2 Nucleocapsid
cross reactivity:	Does not cross react to <i>Influenza</i> type A, <i>Influenza</i> type B, RSV, Adenovirus, RSV, Strep A, Legionella pneumophila, Strep pneumoniae, Mycoplasma pneumoniae, Coronavirus NL63, Coronavirus HKU-1, Coronavirus 229E, Coronavirus OC43
pairing information:	This antibody forms a match pair when used as the capture with A03098 and A03101 as the detection antibodies.
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.