



PRODUCT SPECIFICATION

CATALOG NUMBER: A05010

PRODUCT DESCRIPTION:	Lipocalin-2/NGAL Monoclonal Antibody
IMMUNOGEN:	Recombinant human Lipocalin-2
PURIFICATION:	Protein A or G purified from hybridoma culture supernatant
BUFFER:	PBS, 5% Trehalose, pH 7.4
SPECIFICITY:	Detects human Lipocalin-2.
CROSS-REACTIVITY:	Does not cross-react with recombinant human (rh) Lipocalin-1, rhMMP-9, rmLipocalin-2, or rrLipocalin-2.
STORAGE:	Short term: 2-8°C (1 month) Long term: -20°C (12 months as supplied or 6 months after opening)
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.