



## PRODUCT SPECIFICATION

CATALOG NUMBER: A0501 1

PRODUCT DESCRIPTION:	Lipocalin-2/NGAL Monoclonal Antibody
IMMUNOGEN:	Recombinant human Lipocalin-2/NGAL
PURIFICATION:	Protein A or G purified from hybridoma culture supernatant
BUFFER:	PBS, pH 7.4
SPECIFICITY:	Detects human Lipocalin-2/NGAL. Also recognizes human Lipocalin-2/MMP-9 (Western blot under non-reducing conditions)
CROSS-REACTIVITY:	Does not cross-react with recombinant human (rh) Lipocalin-1 or mLipocalin-2
STORAGE:	Short term: 2-8°C (1 month) Long term: -20°C (12 months as supplied or 6 months after opening)
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.