



## PRODUCT SPECIFICATION

catalog number:	<b>A05016</b>
product description:	NGAL Monoclonal Antibody
lot number:	[lot specific]
immunogen:	Human recombinant NGAL expressed in mammalian cells
source:	Murine
form:	Liquid
purification method:	Protein A chromatography
purity:	≥90% by SDS-PAGE
buffer:	PBS, pH 7.4
preservative:	0.02% NaN <sub>3</sub>
isotype:	IgG <sub>1k</sub>
concentration:	[lot specific] mg/mL
pairing:	This antibody can be used as the capture antibody in conjunction with A05017 (as the detection antibody).
storage:	2-8°C short term (up to three months); -20°C long term (up to 36 months) avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known hazard/biohazard associated with this product.