



PRODUCT SPECIFICATION

catalog number:	A06001
product description:	Renin Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Recombinant human Renin, NSO-derived
source:	Murine
form:	Liquid
purification method:	Protein A or G chromatography; filtered in 0.2 µm solution
buffer:	PBS, pH 7.4
preservative:	none
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	This monoclonal antibody is specific for human Renin, detecting both pro and mature forms in ELISA and Western blot.
cross-reactivity:	<ul style="list-style-type: none">▪ rm Renin 0% ▪ rh Cathepsin D 0%▪ rh Cathepsin E 0%
storage:	Upon receipt store -20°C (up to 12 months); After opening: 2-8°C (one month); -20°C (up to 6 months)
hazard/biohazard:	There is no known hazard or biohazard associated with this product.