



## PRODUCT SPECIFICATION

catalog number:	<b>A06003</b>
product description:	Renin Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human recombinant Renin
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	0.1M PBS, pH 7.4
preservative:	0.5 g/L NaN <sub>3</sub>
isotype:	IgG <sub>1</sub> , Kappa
concentration:	[lot specific] mg/mL
specificity:	Recognizes native and recombinant Human Renin 100% Active Renin                      <10% Prorenin
pairing:	This antibody can be used as the capture antibody in ELISA. Pair with A06004
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.