



## PRODUCT SPECIFICATION

catalog number:	<b>A09167601P</b>
product description:	Elastase Monoclonal Antibody
lot number:	[lot specific]
immunogen:	Elastase CELA3A and CELA3B recombinant Antigens
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	PBS, pH 7.4
preservative:	0.09% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
specificity:	Recognizes the both Elastases
expiration date:	[lot specific]
storage:	-20°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.