



## PRODUCT SPECIFICATION

catalog number:	<b>A09180489P</b>
product description:	Human KL-6 Monoclonal Antibody
lot number:	[lot specific]
immunogen:	Human KL-6
source:	Murine
form:	Liquid
purification method:	Affinity Chromatography
buffer:	15mM Phosphate buffer, 0.15M NaCl
preservative:	0.09% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub> , kappa
concentration:	[lot specific] mg/ml
specificity:	Specific for Human KL-6
pairing:	This antibody can be used in sandwich ELISA as the detection and capture antibody
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.