



PRODUCT SPECIFICATION

catalog number:	A10301
product description:	IL-6 Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Recombinant human IL-6, <i>E. coli</i> derived
source:	Murine
form:	Liquid
purification method:	Protein A/G chromatography; 0.2µm filtered
buffer:	PBS, pH 7.4
preservative:	None
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
endotoxin:	<0.10 EU/µg by LAL
specificity:	This monoclonal antibody is specific for human IL-6.
cross-reactivity:	This antibody does not cross-react with rIL-6, rhOSM, rhLIF, rhIL-11, rhgp130 or rhCNTF.
storage:	Upon receipt store -20°C (up to 12 months); After opening: 2-8°C (one month); -20°C (up to 6 months)
hazard/biohazard:	There is no known hazard or biohazard associated with this product.