



PRODUCT SPECIFICATION

catalog number:	A10308
product description:	IL-6 Monoclonal Antibody
lot number:	[lot specific]
immunogen:	Human Recombinant IL-6
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	PBS 0.1M, pH7.4
preservative:	NaN3 0.5g/L
isotype:	IgG1 kappa
concentration:	[lot specific]
specificity:	Recognizes native and recombinant Human IL-6
pairing:	This antibody can be used in sandwich ELISA as the detector antibody
storage:	2-8°C
expiration date:	[lot specific]
hazard/biohazard:	There is no known biohazard associated with this product.