



## PRODUCT SPECIFICATION

catalog number:	<b>A10380</b>
product description:	Erythropoietin Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	EPO BSA conjugate
source:	Murine
form:	Liquid
purification method:	Protein A chromatography
purity:	>98%
buffer:	PBS
preservative:	0.05% NaN <sub>3</sub>
isotype:	IgG
concentration:	[lot specific] mg/mL
application/notes:	RIA, EIA, immunohistochemistry, and western blot
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.