



## PRODUCT SPECIFICATION

catalog number:	<b>A10462</b>
product description:	AMH (Anti-Mullerian Hormone) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human AMH
source:	Murine ascites
form:	Liquid
purification method:	DEAE chromatography
purity:	>95% by HPLC and SDS-PAGE
buffer:	0.015M potassium phosphate, 0.85% NaCl, pH 7.2
preservative:	0.05% NaN <sub>3</sub>
isotype:	IgG <sub>2a</sub>
concentration:	[lot specific] mg/mL (OD280 nm, E <sup>0.1%</sup> = 1.4)
specificity:	This antibody reacts with human AMH.
pairing:	This antibody can be used in quantitative EIA and paired with A10460503P or A10461503P.
storage:	2-8°C (one year); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.