



PRODUCT SPECIFICATION

catalog number:	A10741
product description:	NT-proBNP (12-16) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Native human proBNP
source:	Murine
form:	Liquid
purity:	>90% pure
buffer:	10mM PBS, pH 7.4
concentration:	[lot specific] mg/mL
specificity:	This monoclonal is specific for the 12-16aa sequence of NT-proBNP.
pairing:	This antibody can be used in ELISA as the conjugate antibody in conjunction with A10751 and/or A10752 (as the capture antibodies).
expiration date:	[lot specific]
storage:	2-8°C short term (two weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.