



PRODUCT SPECIFICATION

catalog number:	A10753
product description:	NT-proBNP (39-76) Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine
form:	Liquid
purity:	>90% pure
buffer:	1XPBS
preservative:	1/1000 P300
concentration:	[lot specific] mg/mL
specificity:	This monoclonal is specific for the 39-76aa sequence of NT-proBNP.
expiration date:	[lot specific]
storage:	2-8°C short term (two weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.