



## PRODUCT SPECIFICATION

catalog number:	<b>A10754</b>
product description:	NT-proBNP (42-46) Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine
form:	Liquid
purity:	>90% pure
buffer:	10 mM PBS, pH 7.4
preservative:	0.1% Proclin300
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
specificity:	This monoclonal is specific for the 42-46aa sequence of NT-proBNP.
expiration date:	[lot specific]
storage:	2-8°C short term (two weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.