



PRODUCT SPECIFICATION

catalog number:	A10760
product description:	NT-proBNP Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Caprylic acid and ammonium purification ; filtered in 0.2 µm solution
buffer:	10mM sodium phosphate, 0.8% NaCl, pH 7.2
preservative:	0.05% NaN ₃
concentration:	[lot specific] mg/mL
specificity:	This monoclonal is specific for the 67-76aa sequence of NT-proBNP.
pairing:	This antibody can be used in ELISA as the capture antibody in conjunction with polyclonal G10800 (as the detection antibody).
manufacture date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.