



PRODUCT SPECIFICATION

| | |
|----------------------|---|
| catalog number: | A10771 |
| product description: | NT-proBNP Monoclonal Antibody |
| lot number: | [lot specific] |
| source: | Murine |
| form: | Liquid |
| purification method: | Protein A chromatography |
| purity: | ≥90% by SDS-PAGE |
| buffer: | PBS, pH 7.4 |
| preservative: | 0.02% NaN ₃ |
| isotype: | IgG _{2b} |
| concentration: | [lot specific] mg/mL |
| titer: | ≥ 1 : 20,000 (ELISA) |
| pairing: | This antibody can be used in sandwich ELISA as the capture antibody in conjunction with A10772 (as the detection antibody). |
| storage: | 2-8°C short term (up to 3 months); -20°C long term (up to 3 years) avoiding multiple freeze thaw cycles. |
| hazard/biohazard: | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |