



PRODUCT SPECIFICATION

catalog number:	A10772
product description:	NT-proBNP Monoclonal Antibody
lot number:	[lot specific]
source:	Murine
form:	Liquid
purification method:	Protein A chromatography
purity:	≥90% by SDS-PAGE
buffer:	PBS, pH 7.4
preservative:	0.02% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
titer:	≥ 1 : 20,000 (ELISA)
pairing:	This antibody can be used in sandwich ELISA as the detection antibody in conjunction with A10771 (as the capture antibody).
storage:	2-8°C short term (up to 3 months); -20°C long term (up to 3 years) avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.