



PRODUCT SPECIFICATION

catalog number:	A11018
product description:	AFP Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human Recombinant AFP
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	0.1 M PBS, pH 7.4
preservative:	0.5 g/L NaN ₃
concentration:	[lot specific] mg/mL
specificity:	Recognizes native and recombinant Human AFP
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.