



PRODUCT SPECIFICATION

CATALOG NUMBER: A11070

PRODUCT:	Monoclonal Antibody to Alpha Fetoprotein (AFP)
IMMUNOGEN:	AFP purified from human fetal cord serum
FORM:	Chromatographically purified in 0.05 M PBS, 0.1% NaN ₃ ; pH 7.4
AFFINITY CONSTANT:	6 X 10 ⁹
SUBCLASS:	IgG
STORAGE CONDITIONS:	2-8°C
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.