



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: A11080**

PRODUCT:	Mouse monoclonal anti-alpha Fetoprotein (AFP)
PURIFICATION:	Protein A
PURITY:	>90% by SDS-PAGE
BUFFER:	10 Mm Phosphate, pH 7.4, containing 150 mM NaCl and 0.1% NaN <sub>3</sub>
SUBCLASS:	IgG <sub>1</sub> , Kappa
STORAGE CONDITIONS:	2-8°C short term -20°C long term