



PRODUCT SPECIFICATION

| | |
|----------------------|---|
| catalog number: | A13021 |
| product description: | <i>Chlamydia trachomatis</i> (LPS) Monoclonal Antibody, Purified |
| lot number: | [lot specific] |
| immunogen: | Chlamydia LGV2 |
| source: | Murine |
| form: | Liquid |
| purity: | >95% by SDS-PAGE |
| buffer: | PBS, pH 7.0 |
| preservative: | 0.02% NaN ₃ |
| isotype: | IgG _{2a} |
| concentration: | [lot specific] mg/mL (based on A ^{1%} _{280nm} = 14.0) |
| storage: | 2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles. |
| hazard/biohazard: | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |