



PRODUCT SPECIFICATION

catalog number:	A13220
product description:	Trichomonas vaginalis Monoclonal Antibody
lot number:	[lot specific]
source:	Murine
form:	Liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	0.01 M PBS, pH 7.2
preservative:	0.1% NaN ₃
isotype:	IgG ₃
concentration:	[lot specific] mg/mL
specificity:	T. vaginalis
pairing:	This antibody can be used in sandwich ELISA as the conjugate antibody in conjunction with A13260041P or A13230041P (as the capture antibodies).
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.