



## PRODUCT SPECIFICATION

catalog number:	<b>A13260041P</b>
product description:	Trichomonas vaginalis Monoclonal Antibody
lot number:	[lot specific]
source:	Murine
form:	Liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	0.01 M PBS, pH 7.2
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG <sub>3</sub>
concentration:	[lot specific] mg/mL
specificity:	T. vaginalis
pairing:	This antibody can be used in sandwich ELISA as the capture antibody in conjunction with A13220041P (as the conjugate antibody).
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.