



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A15010**

|                     |  |
|---------------------|--|
| PRODUCT:            | Monoclonal Antibody to Apolipoprotein A-1  |
| IMMUNOGEN:          | Human Plasma   |
| FORM:               | Immunoglobulin purified from mouse ascites   |
| BUFFER:             | DEAE-chromatography purified IgG in 0.015 M potassium phosphate buffer, 0.15 M NaCl, pH 7.2, with 0.1% NaN <sub>3</sub> as preservative.       |
| SUBCLASS:           | IgG <sub>1</sub>   |
| STORAGE CONDITIONS: | 2-8°C  |
| HAZARD/BIOHAZARD:   | This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product. |