



PRODUCT SPECIFICATION

CATALOG NUMBER: A17030

PRODUCT:	Ferritin Monoclonal Antibody
IMMUNOGEN:	Human Spleen
PURIFICATION:	DEAE
PURITY:	>90% by SDS-PAGE
BUFFER:	10 mM Phosphate, with 150 mM NaCl, pH 7.4, 0.1% NaN ₃
SUBCLASS:	IgG _{2a} , Kappa
STORAGE CONDITIONS:	2-8°C short term -20°C long term
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.