



PRODUCT SPECIFICATION

CATALOG NUMBER: A18080

| | |
|---------------------|---|
| PRODUCT: | hFSH Monoclonal Antibody |
| IMMUNOGEN: | Human FSH, alpha subunit |
| FORM: | Immunoglobulin purified from mouse ascites In 0.05M PBS, pH7.40, with 0.1% NaN ₃ as preservative. |
| AFFINITY CONSTANT: | 4 x 10 ¹⁰ |
| SUBCLASS: | IgG |
| STORAGE CONDITIONS: | 2-8 ⁰ C |
| CROSS REACTIVITY: | Paired with A18090131P. FSH 100%, LH <4%, TSH <1%, Beta-HCG <1%. |
| HAZARD/BIOHAZARD: | This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product. |