



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A19220**

PRODUCT:	HSV 2 Monoclonal Antibody
SOURCE:	Immunoglobulin from mouse ascites
IMMUNOGEN:	Herpes simplex virus (HSV), type 2
FORM:	Each vial contains antibody in PBS, pH 7.4, with 0.1% sodium azide as preservative.
PURIFICATION:	Purified by chromatography on Protein G Sepharose. Purity is tested by electrophoresis.
SUBCLASS:	IgG <sub>1</sub>
STORAGE CONDITIONS:	4°C
REACTIVITY:	Antibody reacts with 56-64 kD band corresponding to Glycoprotein D band in immunoblot of purified HSV.
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.