



PRODUCT SPECIFICATION

CATALOG NUMBER: A19320

ANTIBODY TO:	HSV 1 & 2 Monoclonal Antibody
FORM:	Purified immunoglobulin in 0.2M PB with 0.25M NaCl, pH 7.6, 0.1% sodium azide
SUBCLASS:	IgG _{2a}
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	Highly reactive to HSV-I strains "Mac" and "F". Moderately reactive to HSV-II strain "G". Weakly reactive to HSV-II strain "Ms".
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.