



PRODUCT SPECIFICATION

catalog number:	A20225
product description:	Placental Lactogen (HPL) Monoclonal Antibody
lot number:	[lot specific]
immunogen:	Human Recombinant HPL
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	BBS 0.1M pH 8.0
preservative:	None
isotype:	IgG1
concentration:	[lot specific] mg/mL
specificity:	Recognizes natural and recombinant HPL
pairing:	This antibody can be used in sandwich ELISA as the capture antibody
storage:	2-8 °C
hazard/biohazard:	There is no known biohazard associated with this product.