



PRODUCT SPECIFICATION

CATALOG NUMBER: A21060

PRODUCT:	TSH (TSH & Beta TSH Specific) Monoclonal Antibody
IMMUNOGEN:	TSH
FORM:	Purified immunoglobulin from mouse ascites
BUFFER:	0.05M PBS, 0.1% NaN ₃ , pH 7.4
STORAGE CONDITIONS:	2-8°C
AFFINITY CONSTANT:	2 x 10 ¹⁰
SUBCLASS:	IgG
CROSS REACTIVITY:	TSH 100%, HCG <1%, LH <1%, FSH <1%
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.