



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A21120**

PRODUCT:	TSH, Beta Monoclonal Antibody
IMMUNOGEN:	Human TSH, beta subunit
FORM:	Chromatographically purified Ig fraction in 0.05M PBS, pH 7.4, 0.1% NaN <sub>3</sub>
AFFINITY CONSTANT:	1 x 10 <sup>10</sup>
SUBCLASS:	IgG
STORAGE CONDITIONS:	2- 8°C
CROSS REACTIVITY:	100% TSH, <1% bHCG, Intact HCG <1%, LH <1%, and FSH <1%
COOPERATIVE PAIR:	With A21180
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.