



PRODUCT SPECIFICATION

CATALOG NUMBER: A21180

PRODUCT:	TSH Monoclonal Antibody
IMMUNOGEN:	Human TSH –beta subunit
FORM:	Chromatographically purified Ig fraction in 0.05M PBS, pH 7.4, 0.1% NaN ₃
AFFINITY CONSTANT:	3 x 10 ¹⁰
SUBCLASS:	IgG
CROSS REACTIVITY:	TSH 100% Beta-HCG < 1% Intact-HCG < 1% LH < 1% FSH < 1%
PAIRING INFORMATION:	Paired with A21120
STORAGE CONDITIONS:	2-8°C
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.