



PRODUCT SPECIFICATION

| | |
|----------------------|---|
| catalog number: | A21292 |
| product description: | TSH Monoclonal Antibody, Purified |
| lot number: | [lot specific] |
| immunogen: | Human recombinant TSH |
| source: | Murine |
| form: | Liquid |
| purification method: | Protein A affinity chromatography |
| buffer: | 0.1M PBS, pH 7.4 |
| preservative: | 0.5% NaN ₃ |
| isotype: | IgG ₁ , Kappa |
| concentration: | [lot specific] mg/mL |
| specificity: | Recognizes native and recombinant Human TSH |
| pairing: | This antibody can be used as the detection antibody in ELISA. |
| expiration date: | [lot specific] |
| storage: | 2-8°C |
| hazard/biohazard: | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |