



PRODUCT SPECIFICATION

catalog number:	A21292
product description:	TSH Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human recombinant TSH
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	0.1M PBS, pH 7.4
preservative:	0.5% NaN ₃
isotype:	IgG ₁ , Kappa
concentration:	[lot specific] mg/mL
specificity:	Recognizes native and recombinant Human TSH
pairing:	This antibody can be used as the detection antibody in ELISA.
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.